



# FitNuture Posture Detection App - User Manual

you can often display it in Markdown like this if the Markdown renderer supports it, or you'd add it manually if converting to PDF. For now, we'll assume the app itself displays the logo. -->

## 1. Introduction

Welcome to the FitNuture Posture Detection App! This application is designed to help identify potential postural abnormalities in students using image analysis. By capturing a clear image, the app can provide an initial assessment, generate a report, and store data for future reference and analysis, including potential use in developing more accurate machine learning models.

This manual will guide you through using the app effectively, from capturing the right kind of image to understanding the results and managing data.

## 2. Getting Started

### What You Need:

- A device with a camera (smartphone, tablet, or computer with a webcam).
- Good lighting conditions.
- A clear, uncluttered background.
- The student whose posture is to be analyzed.

### Accessing the App:

Simply navigate to the web address where the FitNuture app is hosted (e.g., your Streamlit Cloud URL).

## 3. Capturing an Accurate Posture Image

The accuracy of the posture analysis heavily depends on the quality of the image you provide. Please follow these guidelines carefully:

### 3.1. Importance of a Good Image

A clear, well-positioned image allows the app's AI to accurately identify key body landmarks, which are crucial for calculating posture metrics. Poor images can lead to inaccurate results or failure to detect a person.

### 3.2. Camera Placement and Subject Positioning

#### a. Distance from Subject:

- Position the camera far enough away so that the entire body of the student, from head to feet, is clearly visible within the frame.
- There should be some space above the head and below the feet.

#### b. Camera Angle and Height:

- The camera should be placed at approximately the student's mid-torso (chest/stomach) height.
- Avoid angling the camera upwards or downwards. It should be as level as possible, pointing directly at the student.

#### c. Full Body Visibility:

- Crucial: The entire body must be in the shot. This includes the top of the head, both arms (even if by the sides), the torso, both legs, and the feet.
- Ensure no body parts are cut off by the edge of the frame.

#### d. Subject's Posture for the Photo:

- The student should stand straight and tall, in a relaxed but upright natural stance.
- They should be facing directly sideways to the camera for a profile view. The app is designed to analyze a sagittal plane (side view) posture. (If a frontal view is required for specific analyses in the future, that would need different instructions).
- Feet should be shoulder-width apart, or a comfortable natural stance, pointing forward.
- Arms should hang naturally at their sides.
- Ask the student to look straight ahead, not up, down, or at the camera (unless the camera is directly in front of them at eye level for a frontal shot, which is not the primary mode for this analysis).

#### e. Clothing:

- Wear form-fitting clothing if possible. Bulky or loose clothing can obscure body landmarks and affect accuracy.
- Avoid clothing with busy patterns that might confuse the landmark detection. Solid, contrasting colors (relative to the background) are best.

#### f. Background:

- Use a plain, uncluttered background. A clear wall is ideal.
- Ensure there are no objects or other people in the background that could interfere with person detection.
- The color of the background should contrast with the student's clothing.

#### g. Lighting:

- Ensure the student is well and evenly lit.
- Avoid shadows falling on the student.
- Natural daylight is often best, but good indoor lighting can also work. Avoid direct, harsh light that creates strong shadows or overexposure.
- Do not have a bright light source (like a window) directly behind the student, as this will make them appear as a silhouette.

### 3.3. Using the App's Camera Input

1. Once on the app page, select the "Use Camera" input mode. 2. If prompted by your browser, allow the app to access your camera. 3. Position the student according to the guidelines above. 4. Use the camera preview in the app to frame the student correctly. 5. Click the "Take a picture" button.

### 3.4. Uploading an Existing Image

If you have a pre-existing image that meets all the guidelines: 1. Select the "Upload Image" input mode. 2. Click the "Browse files" button (or drag and drop an image). 3. Select the image file (JPG, PNG, JPEG) from your device.

## 4. Performing the Analysis

### 4.1. Entering Child's Name

- Before taking a picture or uploading an image, you must enter the child's name in the "Child's Name" field. This is a mandatory field for associating the analysis with the correct individual.

### 4.2. Selecting Abnormalities to Detect

- Below the name input, you'll find a section titled "Select Abnormalities to Detect."
- By default, all listed abnormalities are selected for detection.
- You can uncheck any specific abnormalities you do not want the app to analyze for this particular session.
- There's a "Select All" checkbox to quickly toggle all options on or off.

### 4.3. Understanding the Results

- After a successful image capture and processing, the app will display:
- An image with detected body landmarks overlaid.
- A list of the selected abnormalities and whether they are considered "Present" or "Not Present" based on the app's analysis.
- Important: The results are based on an automated analysis and are for informational purposes only. They are not a substitute for professional medical advice.

## 5. Saving and Managing Results

### 5.1. Saving Results Locally

- After an analysis is complete, if you are satisfied with the result, you can save it.
- Click the "(Save Icon) Save Result Locally" button.
- This will store the analysis data (student details, detected abnormalities, and calculated metrics) within your current browser session. This data is not yet uploaded to the cloud.

### 5.2. Viewing Collected Records

- Scroll down to the "(Chart Icon) View Locally Saved Records" section.
- Here, you'll see a table of all the records you've saved locally during your current session.
- If there are many records, you can use the "Select Page" dropdown to navigate through them.
- You can also use the "Search by Student Name or ID" field to find specific records.

### 5.3. Downloading Local Records (CSV)

- In the "View Locally Saved Records" section, you'll find a "(Download Icon) Download All Local Records (CSV)" button.
- Clicking this will download a CSV (Comma Separated Values) file containing all the records currently saved locally. This file can be opened with spreadsheet software like Microsoft Excel, Google Sheets, etc.

## 6. Uploading Data to the Cloud (Azure SQL)

### 6.1. Purpose of Cloud Upload

- Uploading data to the cloud (Azure SQL database) allows for:
- Permanent storage of analysis records beyond your current browser session.
- Centralized data access for authorized personnel.
- Accumulation of a dataset that can be used for further analysis, research, and the long-term goal of training machine learning models to improve detection accuracy.

### 6.2. How to Upload

1. Ensure you have one or more records saved locally (see section 5.1). 2. Scroll down to the "(Cloud Icon) Cloud Storage (Azure SQL)" section. 3. Click the "(Upload Icon) Upload All Saved Local Records to Azure SQL" button. 4. The app will attempt to connect to the database and upload all records currently stored in the "View Locally Saved Records" table. 5. You will see a success message if the upload is complete, or an error message if issues occur.

### 6.3. What Happens to the Data

- When uploaded, the data for each student (including their ID, name, timestamp, detected abnormalities, and posture metrics) is stored in the secure Azure SQL database.
- If a record with the same "Student ID" already exists in the database, the existing record will be updated with the new information. Otherwise, a new record will be created.

## 7. Generating a PDF Report

After a successful analysis: 1. Click the "(Document Icon) Generate PDF Report" button. 2. The app will create a PDF document containing:

- The FitNurture report title and website.
- The company logo.
- Student's Name, ID, and the analysis Timestamp.
- The landmarked image.
- A summary of detected postural issues.
- General recommendations for any detected conditions.
- A disclaimer.

3. Once generated, a "(Download Icon) Download Report PDF" button will appear. Click this to save the PDF file to your device.

## 8. Troubleshooting

### Common Issues & Solutions:

- "No person detected" or Poor Landmark Detection:
  - Cause: Image quality is likely the issue.
  - Solution: Review section 3.2 ("Camera Placement and Subject Positioning") carefully. Ensure:
    - Full body is visible (head to feet).
    - Good, even lighting (no strong shadows or backlighting).
    - Plain, uncluttered background.
    - Subject is standing straight and sideways to the camera.
    - Appropriate clothing (not too loose).
    - Camera is stable and image is not blurry.
    - Try taking the picture again, adjusting distance or angle slightly.
- PDF Generation Issues (e.g., Blank PDF, Errors):
  - Cause: Could be due to issues with the data being processed or temporary app/browser glitches.
  - Solution:
    - Ensure a successful analysis was completed and results are displayed on the screen.
    - Try refreshing the app page and performing the analysis again.
    - If the problem persists, note any error messages displayed and contact support.
- Cloud Upload Errors (e.g., "Database Connection Error"):
  - Cause: May be due to internet connectivity issues, temporary server problems, or incorrect app configuration (usually an admin concern).
  - Solution:
    - Check your internet connection.
    - Try again after a few minutes.
    - If the issue persists, it might require attention from the app administrator to check database credentials and server status.
- Camera Not Working (especially on Mobile):
  - Cause: Browser might not have permission to access the camera.
  - Solution:
    - When prompted by your browser, ensure you "Allow" camera access for the app's website.
    - Check your browser's settings for camera permissions for the site.
    - Try a different browser if the issue continues.
    - Ensure no other app is currently using the camera.

## 9. Contact / Support

For any issues not covered in this manual or for further assistance, please contact: [www.futurenurture.in](http://www.futurenurture.in) | [info@futurenurture.in](mailto:info@futurenurture.in)

## 10. Disclaimer

This automated analysis is for informational purposes only and not a substitute for professional medical advice. Consult a healthcare provider for health concerns.

\*Thank you for using the FitNurture Posture Detection App!\*